

## Technical Data

Product: DIVAC 1.4 HV3C  
Part No.: 13520V

Inlet connection:	Hose nozzle ID 10, G 1/8" Inside thread	
Outlet connection:	Hose nozzle ID 10, G 1/8" Inside thread	
Pumping speed:	1.3 m <sup>3</sup> /h	0.77 cfm
Ultimate pressure:	≤ 2.0 mbar	≤ 1.5 Torr
Max. outlet back pressure: (absolute)	1500 mbar	1125 Torr
Noise level:	48 dB(A)	
Mains connection:	1-ph, 90-230 V, 50/60 Hz	
Motor rating:	135 W	
Motor rating at ultimate pressure:	35 W	
Rotation speed:	700 - 1600 min <sup>-1</sup>	700 - 1600 rpm
Motor protection rating:	IP 20	
Admissible gas inlet temperature:	+5 to +40°C	+41 to +104°F
Admissible ambient temperature:	+5 to +40°C	+41 to +104°F
Dimensions:	see dimension sheet	
Weight:	8.6 kg	19.0 lbs

Technical data are subject to change