

## Technical Data

Product: DIVAC 2.2 L  
Part No.: 13512

Inlet connection:	Hose nozzle ID 10, G 1/4" Inside thread	
Outlet connection:	Hose nozzle ID 10, G 1/4" Inside thread	
Pumping speed:	2.2 m <sup>3</sup> /h	
Ultimate pressure:	≤ 8 mbar	
Max. outlet back pressure: (absolute)	2000 mbar	
Noise level:	52 dB(A)	
Mains connection:	1-ph, 230 V, ±10%, 50 Hz	
Motor rating:	245 W	
Rotation speed at 50 Hz:	1500 min <sup>-1</sup>	
Motor protection rating:	IP 44	
Admissible gas inlet temperature:	+5 to +40°C	+41 to +104°F
Admissible ambient temperature:	+5 to +40°C	+41 to +104°F
Dimensions:	see dimension sheet	
Weight:	12.6 kg	27.8 lbs

Technical data are subject to change